

# PHONICARE<sup>®</sup>

cleansing wipes



with lemon freshness

clean, hygienic & fresh  
hearing devices



by JSC International



For the surficial care of hearing aids, hearing protection and In Ear Monitoring.

- Removes debris & earwax thoroughly.
- Good maintenance keeps your devices in good shape and provide a longer life.

USAGE: Wipe of debris & ear wax thoroughly but careful and let your device dry.

WARNING: Keep out of the reach of children! In case of contact with eyes: Rinse cautiously with water for several minutes. Do not hold against open flame or other ignition source.

Please read the identification label before use.

CONTAINS: <10% isopropanol CAS 67-63-0), <0.2% anionic surfactants, <0.1% perfume oils

Zur äußeren Pflege von Hörgeräten, Gehörschutz und In Ear Monitoring.

- Entfernt Verschmutzungen & Cerumen gründlich.
- Gute Pflege maximiert die Lebensdauer Ihrer Geräte.

ANWENDUNG: Das zu reinigende Gerät vorsichtig mit dem Tuch abreiben und trocknen lassen. HINWEIS: Darf nicht in die Hände von Kindern gelangen! Bei Berührung mit den Augen: Einige Minuten lang behutsam mit Wasser ausspülen. Nicht gegen offene Flamme oder andere Zündquelle halten. Vor Gebrauch Kennzeichnungsetikett lesen.

ENTHÄLT: <10% Isopropanol (CAS 67-63-0), <0.2% anionische Tenside, <0.1% Parfümöle

JSC International  
Janstuytstraat 82  
2662 JB Bergschenhoek  
The Netherlands  
+31 (0)10 888 3015  
info@phonicare.com

LOT



MD



REF

1016



ECS GmbH  
Wolfener Strasse 32-34  
D-12681 Berlin  
Germany

